Operation Manual

Excelsior Syringe Pump ESP-60

2009 Edition Revision 3

Manual de funcionamiento

Bomba de jeringa ESP-60-ES

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Bomba de seringa ESP-60-PT

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Table of contents

Introduction2

Infusion Methods2

Diagram of the ESP2

Instructions for use.....3

Visual / audible indicators and alarms3

Occlusion correction procedure3

Infusion times and flow rates4

Cleaning4

Periodic functional inspection.....4

Warnings and cautions5

Technical Service

Biomedical Engineering Initial Check-In5

Trouble shooting guide6

Disassembly and replacement of components6

Specifications8

Limited warranty9

Introduction

Excelsior Medical Corp. is a company dedicated to the intermittent administration of drugs via syringe pumps.


The ESP Syringe Pump is a small, compact infusion pump, which infuses doses with 3cc-60cc syringes. It has a plastic case, syringe clamp, pusher block and membrane type keypad. The pump is user friendly - all user information is displayed on the pump's front. Three speeds provide versatility to extend delivery time by a factor of 2 or 3. Only one button to start infusion, no drip counting, activation, spiking or admixture. The pump provides visual and audible alarms for end of infusion, occlusion, low battery, accuracy and system malfunction.

In addition to its syringe pump Excelsior Medical Corp. offers a complete system:

1. Expert analysis of your facility's drug requirements, current method costs and expected cost savings by converting to syringe pumps.
2. State of the art high speed dispensing pump, tubing and adaptors specifically for syringe batching.
3. Experienced personnel for implementation and training to insure realization of the anticipated cost savings without any hassle or headache.
4. Toll-free telephone number: (888) ESP-Pump, (888) 377-7867) for ordering products or any assistance.
5. Highest quality alkaline batteries.

Infusion methods

1. As a primary intermittent infusion system in combination with a heparin lock.
2. As a piggyback in combination with a primary IV delivering compatible maintenance fluids or drug diluents, the primary line flow rate need not be interrupted. The microbore extension set is attached to the primary line at the lower injection port or flashball site.
3. As a piggyback in combination with an IV controller or IV pump. No change to the pump or controller is required. The microbore extension set is attached to the primary line at an injection port or flashball site downstream from the primary infusion device.
4. As a piggyback in combination with a volume control chamber (burette). The flow rate for the burette is not interrupted. The microbore extension set is attached to the primary line at either the lower injection port or the flashball site.

Instructions for use

CDC Guidelines, Institutional Policies and accepted IV therapy principles should be followed.

1. Load batteries. Open the battery compartment by turning the battery screw cap on the bottom of the device. Insert four "C" cell alkaline batteries into the unit with the negative (-) side or flat side of the battery going in first toward the top of the device. *Use only alkaline batteries.*
2. Verify the drug, concentration and dose in the filled syringe are correct for the patient.
3. Using standard aseptic technique and following all guidelines of the hospital, attach the microbore extension set to the filled syringe. Prime the syringe and microbore extension set. *Expel all air making sure to eliminate all air bubbles.*
4. Move the pusher arm mechanism to the top of the device by squeezing the two black levers of the pusher arm together and sliding the pusher arm upward.
5. Place the filled syringe into the syringe barrel clamp holder. *The syringe barrel flange must be fully inserted in the molded slot of the syringe barrel clamp holder.*
6. Move the pusher arm mechanism down until it touches the syringe plunger. *The syringe plunger should be flush against the pusher arm mechanism.*
7. Using standard aseptic technique and following all guidelines of the hospital, attach the distal end of the primed microbore extension set to the infusion site.
8. Start the infusion by pressing the desired *START/SPEED* Key.
"Start Normal "Speed Key-the pusher arm mechanism travels at 5.48 in/hr.
"Start Slow "Speed Key- the pusher arm mechanism travels at 2.74 in/hr.
"Start XSlow "Speed Key-the pusher arm mechanism travels at 1.83 in/hr.
9. The Audio Alarm will automatically be activated (indicated by the green flashing LED on the alarm key) and will sound at the end of infusion or at occlusion. Press the *ALARM* key to disable the Audio Alarm. The green LED will be off if the Audio Alarm is disabled. *The visual end of syringe/occlusion alarm (red flashing LED) cannot be deactivated.*
10. Red flashing LED will indicate the end of infusion. If the Audio Alarm was not disabled, a periodic beep tone will sound. Press the *OFF* key to turn off the alarm and power to the device.
11. To stop the infusion at any time, press the *OFF* key.
12. The device may be attached to an IV pole with the pole clamp or on an IV pole hook by using the stainless steel hanger on the top of the unit. The pole clamp is firmly attached to the unit by a tightened set screw on the side of the pole clamp aluminum mounting bracket.

Visual /audible indicators and alarms

1. Flashing green LED above one of the three *START/SPEED* keys indicates the pump is operating at that selected speed. To change the selected speed of operation,
 - a) press the *OFF* key to stop infusion
 - b) press the desired *START/SPEED* key to begin infusion
2. Flashing green LED on the *Audio Alarm* key indicates that the *Audio Alarm* is activated and will sound a periodic beep tone at the end of infusion or occlusion. To disable the Audio Alarm, press the *Audio Alarm* key. The flashing green LED on the *Audio Alarm* key will be off.
3. Flashing red LED (next to "Attention" on the keypad) indicates an end of infusion or occlusion. If the syringe plunger is at the bottom of the syringe, the alarm indicates the end of infusion. Otherwise, the alarm indicates an occlusion. It also flashes when the low battery light is activated.
4. Flashing orange LED (next to "Low Battery") and flashing red LED indicate the batteries have low power. If the low battery light appears during an infusion, the infusion may safely be continued to completion. Prior to starting a new infusion, replace all four batteries with new "C" cell alkaline batteries. The ESP will not start a new infusion if the low battery LED light is flashing.
5. All Speed lights and Attention light flashing indicate unit accuracy has exceeded +/- 3% error. Check to see if the pusher arm was prevented or restricted from movement. Correct the pusher arm obstruction and continue the infusion or contact biomedical engineering if the alarm continues.

Occlusion correction procedure

1. Press the *OFF* key to turn off the *Audio Alarm* and flashing red LED.
2. Relieve the pressure by moving the Pusher arm toward the top of the unit away from the syringe plunger and/or Using standard aseptic technique and following all guidelines of the hospital, disconnect the microbore extension set at the infusion site to prevent an accidental bolus.
3. After clearing the occlusion, move the pusher arm mechanism down until it touches the syringe plunger. *The syringe plunger should be flush against the pusher arm mechanism.*
4. Using standard aseptic technique and following all guidelines of the hospital, attach the distal end of the primed microbore extension set to the infusion site.
5. Start the infusion by pressing the desired *START/SPEED* key. The *Audio Alarm* will automatically be activated

Infusion times and flow rates

(Monoject® Syringes)

Syringe		Infusion Times (in minutes)			Flows Rates mL/hr			Nominal Occlusion
Volume	Size	X-S	Slow	Normal	X-S	Slow	Normal	PSI
5mL	6mL	51	34	17	5.9	8.8	17.7	41
10mL	12mL	66	44	22	9.2	13.7	27.4	26
20mL	20mL	80	53	27	15.1	22.6	45.2	16
30mL	35mL	87	58	29	20.7	31.0	62.0	12
50mL	60mL	116	78	39	25.8	38.6	77.2	9.3

(Becton Dickinson Syringes)

Syringe		Infusion Times (in minutes)			Flows Rates mL/hr			Nominal Occlusion
Volume	Size	X-S	Slow	Normal	X-S	Slow	Normal	PSI
5mL	5mL	57	38	19	5.2	7.8	15.7	45
10mL	10mL	78	52	26	7.6	11.4	22.7	31
20mL	20mL	90	60	30	13.3	19.8	39.7	18
30mL	30mL	105	70	35	17.1	25.6	51.1	14
50mL	60mL	116	77	39	25.9	38.8	77.5	9.2

Cleaning

Clean the exterior surfaces of the unit with a damp cloth and mild detergent or isopropyl alcohol. Avoid getting liquids inside the device.

CAUTION: Use of harsh cleaning agents or solvents may cause damage to the plastic case and components.

Periodic functional inspection

A periodic functional check should be performed at least every six months. The inspection should include the following:

1. Syringe barrel clamp holder should adequately hold each size of syringe (5/6, 10/12, 20, 30/35, 60) firmly in place.
2. The pusher arm mechanism should release by pressing the two black release levers and move easily up and down the unit.
3. The pusher arm mechanism should engage after releasing the two black release levers.
4. Press the START NORMAL speed key. The green LED directly above that key should be flashing and motor operation should be slightly audible.
5. Place the pusher arm mechanism completely at the bottom of travel next to the syringe barrel clamp holder. The operation should continue until the end of infusion alarm. Both the flashing red LED and the periodic beep tone should be operative.
6. Press the ALARM key. The audible periodic beep tone should cease, but the flashing red LED should continue.
7. Press the OFF key and then the START NORMAL speed key. The unit should not restart, but continue to indicate both visual and audible end of syringe alarms.
8. Repeat steps 4 and 5 for the two other speeds.
9. Pusher arm mechanism and syringe barrel clamp should move freely on the two rails. Use petroleum jelly to lightly lubricate the two rails and lead screw. **DO NOT USE SILICONE LUBRICANTS - DAMAGE TO PLASTIC MAY OCCUR.**

Warnings and cautions

- Avoid any liquids inside the pump. The unit is not waterproof or water resistant.
- If a large amount of fluid is spilled on or in the unit, send to the hospital biomedical department prior to performing an infusion.
- Do not sterilize the pump.
- Do not use the pump with narcotics, anesthesia, chemotherapy, radioactive pharmaceuticals, fat emulsions or blood products.
- This unit is intended for use under hospital supervision only. Home patient use is allowed only with instructions provided to the home patient.
- Use only with supervision and direction from a licensed pharmacist, nurse or physician.
- Do not reuse disposable syringes or microbore extension sets. Both are for Single Use Only. Follow CDC and or institutional guidelines. Use only tubing with Excelsior specifications.
- Battery life will be affected by alarm duration and backpressure.

Technical Service

Biomedical Engineering Initial Check-in Procedure for the Excelsior Syringe Pump

Batteries are already installed at the factory prior to shipment. The unit uses 4 "C" cell alkaline batteries with the (-) flat side toward the top or into the unit. (+) side is toward the battery cap.

Syringe clamp holder should adequately hold each size of syringe (5/6, 10/12, 20, 30/35, 60) firmly in place.

Pusher arm mechanism should release by pressing the two black levers and move easily up and down the Unit. The pusher arm mechanism should engage after releasing the black levers.

Visual and audible indicators/alerts

1. Press the START NORMAL SPEED key. The green LED directly above that key should be flashing and motor operation should be slightly audible.
2. Place the pusher arm mechanism completely at the bottom of travel next to the syringe barrel clamp holder. The operation should continue until the end of infusion alarm. Both the flashing red LED and the periodic beep tone should be operative.
3. Press the ALARM key. The audible periodic beep tone should cease, but the flashing red LED should continue.
4. Press the OFF key and then the START NORMAL SPEED key. The unit should not restart, but continue to indicate both visual and audible end of syringe alarms. To restart the pump the pressure must be relieved by moving the pusher block arm from the syringe.
5. Repeat steps 4 and 5 for the two other speeds.

Linear Accuracy is verified up to 670 times per second by an optical encoder. If the motor is not operating within specifications, all three speed lights and the attention LED will flash simultaneously. To test the accuracy of the pump, you need to measure the rate of travel of the pusher block mechanism. Measurement of other parameters such as volume will add additional variables (and potential error) into the test. The following two procedures are recommended.

1. Insert a filled syringe into a pump and compare the infusion time to that indicated in the "Infusion times and flow rates" chart in the ESP user manual.
Note: Do not include pusher block travel beyond the end of syringe in any test for accuracy of the pump. The time it takes the pusher block mechanism to travel from the end of the syringe to the trip point that activates the attention alarm varies from pump to pump.
2. Insert a calibrated dial caliper between the flat surface of the pusher block mechanism and the flat surface at the end of the syringe plunger. Measure the distance traveled by the pusher block over a period of time. The linear rate of travel should be equivalent to:
 - a) 5.48 inches per hour +/- 3% at Normal Speed
 - b) 2.74 inches per hour +/- 3% at Slow Speed
 - c) 1.83 inches per hour +/- 3% at XSlow Speed

Occlusion pressure can be tested with the following procedure:

1. Fill a 20ml Monoject syringe with 20ml of water. Connect the micro-bore tubing to the syringe and prime the line.
2. Connect the micro-bore to a pressure meter and purge any air in the system (use of a three-way stopcock may be helpful).
3. Insert the syringe into the ESP. Push the desired Start/Speed key to start the infusion.
4. When the ESP detects an occlusion, the red "Attention" LED will flash simultaneously with the audible alarm and the motor will stop. Pressure readings should be in the range of 8psi - 24psi with a nominal value of 16psi.

Note:

- a) A new syringe should be used for each test.
- b) If another syringe size is used for the test, the tolerance from the nominal pressure PSI listed in the "Infusion times and flow rates" chart is plus or minus fifty percent ($\pm 50\%$).
- c) PSI is determined by dividing the constant occlusion force of 8.0 +/- 1.0 pounds by the cross sectional area (in square inches) of each syringe. Therefore, a larger cross sectional area (i.e. 60mL syringe) yields a lower nominal occlusion psi than a smaller cross sectional area (i.e. 12mL syringe).

Trouble Shooting Guide

No Power to Unit

- No batteries or dead batteries
- Improper battery installation (batteries must be installed with "+" toward bottom)
- Battery screw cap is not installed
- Metal contact on battery screw cap is missing
- Poor contact between battery cap and tab
- Disconnected/damaged wire from battery compression spring ("- end) to circuit card.
- Disconnected/damaged wire from battery cap contact spring ("+" end) to circuit card
- Incorrect connection of keypad to circuit card
- Inoperative membrane keypad.
- Malfunctioning circuit board.

No infusion - LED's flashing

- Pusher arm compression spring is missing
- Lead screw is damaged or disconnected from motor or lead screw mounts
- Pusher arm and/or half-nuts are damaged
- Damaged or loose lead screw gear. Missing or not fully seated E-clip on lead-screw. Lead-screw is not properly lubricated.
- Damaged motor or motor gear-head
- Motor is not properly attached to motor mount.

LED's not flashing and/or membrane keys not operative

- Disconnected membrane keypad
- Damaged LED's/membrane keypad
- Malfunctioning circuit board

No Audio and Visual Alarms at end of infusion or occlusion

- Occlusion slide potentiometer is not connected to syringe clamp slider. The slider is not seated properly (fully) on the occlusion slide potentiometer
- Damaged or loose gear or lead- screw. Missing or not fully seated E-clip on lead-screw. Lead-screw is not properly lubricated.
- Damaged pusher arm, half-nuts or lead-screw
- Improper wire connection from occlusion slide potentiometer to circuit card
- Circuit board malfunction.

No Audio Alarm - Visual Alarm operates

- "Alarm" key was used to turn off audio alarm.
- Damaged buzzer

Pusher arm is difficult to move on rails

- Damaged guide rails or lead-screw
- Damaged pusher arm.
- Rails not properly lubricated

Disassembly and Replacement of Components

Opening the unit

1. Remove the battery screw cap and batteries.
2. Remove the stainless steel wire hangar from the top of the unit by expanding each side outward.
3. Place the unit on a flat surface with the keypad facing downward.
4. Use a Phillips screwdriver to remove the two flat-head 6/32 x 1 3/4" machine screws that hold the pole clamp mount to the case.
5. Use a Phillips screwdriver to remove the remaining two flat-head 6/32 x 1" machine screws that hold the halves of the case together.

CAUTION: Keep the longer screws with the pole clamp mounting bracket. Installing screws that are longer than one inch without the pole clamp mount will cause damage to the keypad.

Closing the unit

1. Place the front (keypad) of the unit on a flat surface with the keypad facing downward.
- CAUTION: The surface must be free of debris to prevent damage to the keypad.**

2. Connect the wire from the battery compression spring wire ("_" end) to the circuit card at pin location J3.
3. Connect the wire from the battery contact spring ("+" end) to the circuit card at pin location J2.
4. Connect the keypad ribbon cable to the circuit card at pin location J5.
5. Install the mechanical lead- screw/guide rail subassembly into the case. The motor mount and rear mount should easily install onto the two sets of plastic guide pins in the case.
6. Connect the motor cable to the circuit card at pin location J4. Ribbon cable points away from buzzer.
7. Make sure the motor cable is out of way; do not close the back cover on the motor cable.
8. Check all wire connections are proper and secure.
9. Ensure that the syringe clamp slider is fully grabbing the occlusion slide potentiometer the black moveable slide piece.

CAUTION: Black slide of occlusion slide potentiometer must be securely engaged with the syringe clamp slider mechanism for proper occlusion detection.

Replacement of syringe clamp holder

1. Use a Phillips screwdriver to remove the two flat-head 8/32 x 3/8" machine screws that hold the base of the syringe clamp to the slide mechanism.
2. Position the new syringe clamp holder on the slide mechanism and install a screw in the hole toward the top of the unit.
3. Hold the top of the syringe barrel clamp open to allow access to the second hole. Install the second screw .

Replacement of lead- screw rail assembly

1. Disconnect the motor cable from the circuit card at location J4,
2. Remove the lead screw rail assembly from the front case
3. Remove the motor from the lead-screw rail assembly and put aside. See below replacement of motor.
4. Lightly lubricate the new lead-screw and rails with petroleum jelly.

CAUTION: Do not use silicone lubricants. Damage to the plastic components may occur.

5. Check that the gear or the lead screw rotates smoothly with slight force from the hand.
6. Verify E-clip is in place and/or secure
7. Install the motor onto the new lead-screw rail assembly. Do not over tighten the screws.
8. Install the new lead-screw rail assembly with motor onto the front case and connect the motor cable to the circuit card at location J4. The stripe of the motor cable is facing down.

Replacement of motor

CAUTION: Do not turn the motor gear by hand

1. Disconnect motor cable from circuit card at location J4.
2. Remove the lead screw rail assembly from the front case
3. Use a Phillips screwdriver to remove the three pan-head M2 x 6mm machine screws that hold the motor to the motor mount.
4. Pull the motor with gear through the hole in the motor mount.
5. Check that the gear or the lead screw rotates smoothly with slight force from the hand.
6. If needs lightly lubricate the new lead-screw and rails with petroleum jelly.

CAUTION: Do not use silicone lubricants. Damage to the plastic components may occur.

7. Insert the new motor with gear through the hole in the motor mount. Orient the motor such that the ribbon cable is facing away from the black half nuts. The stripe of the motor tail should be in the down direction.
8. Use a Phillips screwdriver to install the three pan-head M2 x 6mm machine screws that hold the motor to the motor mount. Do not over tighten the screws.
9. Install the new lead-screw rail assembly with motor onto the front case and connect the motor cable to the circuit card at location J4. The stripe of the motor cable is facing down.

Replacement of pole clamp

1. Loose the set screw with a hex wrench.
2. Slide the pole clamp downward off the pole clamp mount.
3. Slide the new pole clamp upward onto the pole clamp mount.
4. Tighten the set screw with a hex wrench only to the point that the pole clamp is immovable.

Replacement of pole clamp Mount

1. Remove the pole clamp from the mount. See above.
2. Use a Phillips screwdriver to remove the two flat-head 6/32 x 1 3/4" machine screws that hold the pole clamp mount to the case.
3. Attach the new pole clamp mount.

Replacement of hanger wire

1. Remove the stainless steel wire hanger from the top of the unit by expanding each side outward.
2. Install the new stainless steel wire hanger by expanding it and inserting each end into the holes on the top of the unit.

Replacement of rear instruction label

1. Peel off the old label. Use a heat gun to loosen the adhesive if necessary.
2. Carefully peel back only 1/2" - 1" of the paper backing on the new label and insert the serial number.
3. Align the label on the back of the unit and press the label section without paper backing firmly onto the case.
4. Slowly pull the paper backing from under the label. Ensure that no air bubbles are formed and the label is properly aligned.
5. Press firmly all around the label for adhesion.

Replacement of occlusion slide potentiometer

1. Disconnect the motor cable.
2. Remove the lead- screw rail assembly from the case.
3. Disconnect the occlusion slide potentiometer wire connector from the circuit card at location J1.
4. Pry the old slide potentiometer off the four tabs of the case. **Be careful not to bend or damage the tabs.**
5. Press the new slide potentiometer firmly onto the four tabs of the case.
6. Connect the wire connector to the circuit card at location J1.

Note: From the view of looking at the component side of the circuit card, the wire orientation must be with the green wire from the left side of the slide potentiometer matching up with left pin, the right with the right and middle with the middle of the circuit card at location J1.

Specifications

Accuracy (Travel): +/- 3%

Battery Life: 400 hours typical

Low Battery Level: If low battery LED is flashing - replace all four batteries with new alkaline "C" cells

Batteries: Four "C" size alkaline cells

Syringes Recommended: Monoject single use 6ml - 60ml, BD or Terumo
(3mL for saline flushes)

Rate of Travel:

Normal 5.48 +/- 3% inches per hour

Slow 2.74 +/- 3% inches per hour

XSlow 1.83 +/- 3% inches per hour

Occlusion Force: 8.0 +/- 1.0 pounds (refer to chart for syringe psi; syringe psi is determined by occlusion force and syringe diameter)

Overall size: 9.2" x 1.7" x 4.7" (approximate)

Weight: 1.2 lbs not including batteries

UL Classified Symbols on

Unit:

Refer to documentation for additional information (refer to information on this page)

Type B protection against shock.

Classified by Underwriters Laboratories Inc. for Canada

Classified by Underwriters Laboratories Inc. for USA

EXCELSIOR SYRINGE PUMP
CLASSIFIED BY UNDERWRITERS LABORATORIES INC.
WITH RESPECT TO ELECTRICAL SHOCK, FIRE, MECHANICAL
AND OTHER SPECIFIED HAZARDS ONLY IN ACCORDANCE WITH
UL2601-1 AND CAN/CSA

Limited warranty

Product Standards: Excelsior products are certified to comply with written Excelsior specifications at the time of shipment.

Product Warranty: The Excelsior Syringe Pumps are warranted against defects in material and workmanship for a period of six months from the date of delivery. Excelsior will repair or replace, at its option, any product that proves to be defective during the warranty period, provided proper use and maintenance procedures have been followed as described in the operation manual. If upon examination, Excelsior determines that misuse of product is the cause for the repair, all labor, material and shipping costs involved shall be paid by the Buyer.

All defective products or components shall be returned to Excelsior with a detailed explanation of the failure. Buyer shall obtain a return authorization prior to return of product and all transportation charges must be prepaid.

Limitations: THE PROVISIONS ABOVE CONSTITUTE EXCELSIOR'S SOLE OBLIGATION AND EXCLUDE ALL OTHER REMEDIES OR WARRANTIES EXPRESS, STATUTORY OR IMPLIED, INCLUDING WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, whether or not purposes or specifications are described herein. Excelsior further disclaims any responsibility whatsoever to Buyer or to any other person for injury to person or damage to or loss of property or value caused by any product that been subjected to misuse; negligence or accident; or misapplied or used in violation of product manuals, instructions or warnings; or modified or repaired by unauthorized persons; or improperly installed.

UNDER NO CIRCUMSTANCES SHALL EXCELSIOR BE LIABLE FOR ANY INCIDENTAL, INDIRECT, CONSEQUENTIAL OR SPECIAL DAMAGES, LOSSES OR EXPENSES (even if Excelsior had been advised of the possibility of such damages) arising from the use of, inability to use or repair of Excelsior's product for any purpose whatsoever. BUYER AGREES THAT ANY RECOVERY AGAINST EXCELSIOR SHALL NOT BE GREATER THAN THE PURCHASE PRICE PAID FOR THE PRODUCT AS TO WHICH SUCH CLAIM IS MADE. EXCELSIOR'S SELECTION TO ONE OF ITS ALTERNATIVES (REPLACE OR REPAIR) SHALL BE BUYER'S EXCLUSIVE AND SOLE REMEDY. No person is authorized to modify or amend the written warranty, product claims and specifications or to bind the corporation to any term contrary to the terms herein.

Service performed by other than Excelsior or its authorized agents may, at the discretion of Excelsior, be cause to void this warranty. Excelsior's tubing sets are the only tubing products that may be used with the ESP pump. If other sets are used, it may, at the discretion of Excelsior, be cause to void this warranty.

Excelsior Medical Corp
1933 Heck Ave
Neptune, NJ 07753 USA
(800) 4-U-PHARM
www.excelsiormedical.com

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Classificadas pelo Underwriters Laboratories Inc. para o Canadá
Classificadas pelo Underwriters Laboratories Inc. para E.U.A.